


Neopuff™ INFANT RESUSCITATOR

PRODUCT CATALOGUE



	Section
Neopuff™ Infant Resuscitator	1
Neopuff™ Infant Resuscitator (Stand Alone)	1.1
Neopuff™ Infant Resuscitator (Infant Warmer Integrated Module)	1.2
Neopuff™ T-Piece and Resuscitation Circuit	2
Single Use Adjustable T-Piece	2.1
Humidified Resuscitation Circuit	2.2
Reusable T-Piece	2.3
Reusable Patient Supply Line	2.4
Neonatal Resuscitation Masks	3
Single Use Neonatal Resuscitation Masks	3.1
Resuscitation Kits	4
Single Use Resuscitation Kits	4.1
Reusable Resuscitation Kits	4.2
Neopuff™ Accessories	5
Reusable Test Lung	5.1
Reusable Gas Supply Line	5.2
Gas Blenders	5.3
Humidifier	5.4
Neopuff™ Mounting Options	6
Pole Mount	6.1
Wall Mount	6.2
Pole and Rail Central Mount	6.3
Part Number Index	7

① Neopuff™ Infant Resuscitator Stand Alone




The Neopuff™ Infant Resuscitator is a light weight, stand alone resuscitator unit suitable for Delivery Suite, Post Natal Wards, Community Theaters, Special Care Baby units and Neonatal Intensive Care units. Current resuscitation guidelines recommend the use of T-Piece device with the potential benefits of achieving controlled target PIP and delivering consistent PEEP to help establish FRC and improve lung volume.

The Neopuff™ Infant Resuscitator is easy to use ensuring that irrespective of users' experience, training and fatigue level, the baby will receive safe and consistent breath patterns.

Stand Alone Unit

Product Code	Description
RD900ADU	Neopuff™ Infant Resuscitator - DANISH
RD900AEU	Neopuff™ Infant Resuscitator - ENGLISH
RD900AFU	Neopuff™ Infant Resuscitator - FRENCH
RD900AGU	Perivent™ Infant Resuscitator - GERMAN
RD900AJU	Neopuff™ Infant Resuscitator - JAPANESE
RD900ALU	Neopuff™ Infant Resuscitator - ITALIAN
RD900AMU	Neopuff™ Infant Resuscitator - FINNISH
RD900ANU	Neopuff™ Infant Resuscitator - DUTCH
RD900ASU	Neopuff™ Infant Resuscitator - SPANISH
RD900ATU	Neopuff™ Infant Resuscitator - PORTUGUESE
RD900AVU	Neopuff™ Infant Resuscitator - NORWEGIAN
RD900AWU	Neopuff™ Infant Resuscitator - SWEDISH
RD900AZU	Neopuff™ Infant Resuscitator - CHINESE

② Neopuff™ Infant Resuscitator IW Module




The Neopuff™ Infant Resuscitator is also available as an integrated module inserted into the front of the CosyCot™ Infant Warmer.

An added feature to the CosyCot™ Infant Warmer, this allows easy access for the operator in resuscitating the patient, especially during transport.

Integrated Neopuff™ Module


Product Code	Description
900IW130D	Neopuff™ Infant Resuscitator IW Module - DANISH
900IW130E	Neopuff™ Infant Resuscitator IW Module - ENGLISH
900IW130F	Neopuff™ Infant Resuscitator IW Module - FRENCH
900IW130G	Perivent™ Infant Resuscitator IW Module - GERMAN
900IW130J	Neopuff™ Infant Resuscitator IW Module - JAPANESE
900IW130L	Neopuff™ Infant Resuscitator IW Module - ITALIAN
900IW130M	Neopuff™ Infant Resuscitator IW Module - FINNISH
900IW130N	Neopuff™ Infant Resuscitator IW Module - DUTCH
900IW130S	Neopuff™ Infant Resuscitator IW Module - SPANISH
900IW130T	Neopuff™ Infant Resuscitator IW Module - PORTUGUESE
900IW130V	Neopuff™ Infant Resuscitator IW Module - NORWEGIAN
900IW130W	Neopuff™ Infant Resuscitator IW Module - SWEDISH
900IW130Z	Neopuff™ Infant Resuscitator IW Module - CHINESE

**1 Single Use Adjustable T-Piece**

The adjustable T-Piece is provided with a valve that can be manipulated to achieve desired PEEP level.


The single use adjustable T-Piece includes connected patient supply line.

Product Code	Description	
900RD010	Single Use Adjustable T-Piece	10/bx

**2 Single Use Humidified Resuscitation T-Piece**


The Fisher & Paykel Humidified Resuscitation T-Piece with adjustable PEEP valve is provided with a heated circuit designed to deliver warm, humidified gas to an infant during resuscitation. It is compatible with MR730 and MR850 Fisher & Paykel Humidifiers.

Product Code	Description	
900RD110	Single Use Humidified Resuscitation Kit	10/bx

**3 Reusable Patient T-Piece with PEEP Valve**

The Reusable T-Piece with adjustable PEEP valve is autoclavable. It allows 15mm F connection.

Product Code	Description	
500RD104	Reusable Patient T-Piece with PEEP Valve	


**4 Reusable Patient Supply Line**

The Reusable Patient Supply Line is provided with silicone cuffs and is autoclavable.

Product Code	Description	
500RD103	Reusable Patient Supply Line - 10mm x 1.2m	


1 Single Use Neonatal Resuscitation Masks

The Fisher and Paykel Neonatal Resuscitation Mask is easy to use. It is designed to facilitate an ideal seal on the infants face. The wide size range includes masks for premature neonates.



Product Code	Description	
RD800-EN	Neonatal Resuscitation Mask Starter Kit	5 sizes/bx
RD803-10	Neonatal Resuscitation Mask 35mm (1.38in)	10/bx
RD804-10	Neonatal Resuscitation Mask 42mm (1.65in)	10/bx
RD805-10	Neonatal Resuscitation Mask 50mm (1.97in)	10/bx
RD806-10	Neonatal Resuscitation Mask 60mm (2.36in)	10/bx
RD807-10	Neonatal Resuscitation Mask 72mm (2.83in)	10/bx


① Single Use Resuscitation Kits



The Adjustable T-Piece Resuscitation Kit includes a neonatal resuscitation mask. There are 3 combinations to suit the 3 different mask sizes.

Product Code	Description	
900RD014-10	Adjustable T-Piece + 42mm Mask	10kits/bx
900RD015-10	Adjustable T-Piece + 50mm Mask	10kits/bx
900RD016-10	Adjustable T-Piece + 60mm Mask	10kits/bx

② Reusable Resuscitation Kits



Contains autoclavable patient supply line, Patient T-Piece and gas supply line.


- Lightweight Patient T-Piece for easy operation
- Adjustable PEEP Valve on Patient T-Piece for fingertip control
- Autoclavable

Product Code	Description
500RD002	Reusable Resuscitation Kit
Includes	
500RD103	Reusable Patient Supply Line
500RD104	Reusable Patient T-Piece
500RD107	Reusable Gas Supply Line

**1 Reusable Test Lung**


The reusable test lung is made of silicone and is autoclavable. The test lung is latex free.

Product Code	Description
RD020-01	Reusable Test Lung 1/bx

**2 Reusable Gas Supply Line**


- Length 180mm (7") with 10mm M gas inlet connection

Product Code	Description
900RD008	Gas Supply Line (Short) 10/bx




- Length 1.8m (6') with 10mm M gas inlet connection


Product Code	Description
900RD009	Gas Supply Line (Long) 10/bx



Blender
DISS & NIST
fittings available



Bird MBlender Pole
Mount Bracket



3 Gas Blenders

Gas Blenders are available on an individual basis.


They come complete with color coded nipples and control knobs.

- Flowrate 2-16 L/min (Bird)
0-15 L/min (Bio-Med)
- Includes flowmeter nipples (green or white)


Note: Requires 900IW320 Bird Bracket Blender Complete to be fitted to the CosyCot™ Infant Warmer. Gas hoses not included. Please refer to Fisher & Paykel Infant Warming Catalogue

Product Code	Description
BC501-01	Bird MBlender 2-16LPM DISS
BC502-01	Bird MBlender 2-16LPM NIST
BC503-01	Blender Oxy/Air Bird (AUST)
BC506-01	Blender Oxy/Air Bio-Med (DISS) green
BC507-01	Blender Oxy/Air Bio-Med (DISS) white
BC508-01	Blender Oxy/Air Bio-Med (NIST)
BC510-01	Bird MBlender Pole Mount Bracket

4 Humidifier




MR850 Humidifier




HUMIDIFICATION	
Product Code	Description
MR850	Humidifier
900MR801	Heater Wire Adaptor MR850
900MR823	Heater Wire/Temp Probe Kit
900MR868	Probe Temp/Flow Ang 1450mm

For additional information on humidification devices and accessories, please refer to the Fisher & Paykel Humidification Catalogue or contact your Fisher & Paykel Representative.



1 Pole Mount



The Neopuff™ Infant Resuscitator can be mounted on a pole using a side mounting block and a C-Clamp or a Pole Bracket





Product Code	Description
900RD301	Side Mounting Block
900MR170	C-Clamp
900MR030	Pole Bracket

2 Wall Mount

The Neopuff™ Infant Resuscitator can be mounted on the wall using 2 mounting tongues and a wall mount bracket.




Product Code	Description
900RD302	Wall Mount Bracket (with 2 mounting tongues)

③ Pole and Rail Central Mount


For mounting the Neopuff™ Infant Resuscitator centrally on a pole or on standard rails (2.5-5.5cm x 1cm / 0.98-2.17in x 0.39in), a central mounting block can be affixed to the back of the unit and will fit the Fisher & Paykel bracket and clamps.

Product Code	Description
RD050-01	Neopuff Central Mounting Block Includes: 1 mounting tongue 1 mounting bar 2 M6 x 16 csk screws 4 M6 btn head screws


Product Code	Description
900MR170	C-Clamp
900MR030	Pole Bracket
900MR088	Rail Bracket



C-Clamp



Rail Bracket



Pole Bracket

Part Number	Description	Section
Neopuff™ Infant Resuscitator		1
RD900	Neopuff™ Infant Resuscitator Stand Alone	1.1
900IW130	Neopuff™ Infant Resuscitator Integrated Module	1.2
Neopuff™ T-Piece and Resuscitation Circuit		2
900RD010	Single Use Adjustable T-Piece 10/bx	2.1
900RD110	Single Use Humidified Resuscitation Kit 10/bx	2.2
500RD104	Reusable Patient T-Piece with PEEP Valve	2.3
500RD103	Reusable Patient Supply Line - 10 mm x 1.2m	2.4
Neonatal Resuscitation Masks		3
RD800-EN	Neonatal Resuscitation Mask Starter Kit 5sizes/bx	3.1
RD803-10	Neonatal Resuscitation Mask 35mm (1.38in) 10/bx	3.1
RD804-10	Neonatal Resuscitation Mask 42mm (1.65in) 10/bx	3.1
RD805-10	Neonatal Resuscitation Mask 50mm (1.97in) 10/bx	3.1
RD806-10	Neonatal Resuscitation Mask 60mm (2.36in) 10/bx	3.1
RD807-10	Neonatal Resuscitation Mask 72mm (2.83in) 10/bx	3.1
Resuscitation Kits		4
900RD014-10	Adjustable T-Piece + 42mm Mask 10kits/bx	4.1
900RD015-10	Adjustable T-Piece + 50mm Mask 10kits/bx	4.1
900RD016-10	Adjustable T-Piece + 60mm Mask 10kits/bx	4.1
Neopuff™ Accessories		5
500RD002	Reusable Resuscitation Kit	5.1
RD020-01	Reusable Test Lung 1/bx	5.1
900RD008	Gas Supply Line (Short) 10/bx	5.2
900RD009	Gas Supply Line (Long) 10/bx	5.2

Part Number	Description	Section
BC501-01	Bird MBlender 2-16LPM DISS	5.3
BC502-01	Bird MBlender 2-16LPM NIST	5.3
BC503-01	Blender Oxy/Air Bird (AUST)	5.3
BC506-01	Blender Oxy/Air Bio-Med (DISS)	5.3
BC507-01	Blender Oxy/Air Bio-Med (DISS)	5.3
BC508-01	Blender Oxy/Air Bio-Med (NIST)	5.3
BC510-01	Bird MBlender Pole Mount Bracket	5.3
MR850	Humidifier	5.4
900MR801	Heater Wire Adaptor MR850	5.4
900MR823	Heater Wire/Temp Probe Kit	5.4
900MR868	Probe Temp/Flow Ang 1450mm	5.4
Neopuff™ Mounting Options		6
900RD301	Side Mounting Block	6.1
900MR170	C-Clamp	6.1
900MR030	Pole Bracket	6.1
900RD302	Wall Mount Bracket (with 2 mounting tongues)	6.2
RD050-01	Neopuff Central Mounting Block	6.3
900RD301	Side Mounting Block	6.3
900MR170	C-Clamp	6.3
900MR030	Pole Bracket	6.3
900MR088	Rail Bracket	6.3

INTERNATIONAL P O Box 14 348, Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com
AUSTRALIA Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **AUSTRIA** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **BENELUX** Tel: +31 78 644 0924 Fax: +31 78 644 0914
CHINA Tel: +86 20 3205 3486 Fax: +86 20 3205 3495 **FRANCE** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **GERMANY** Tel: +49 7182 93777 0 Fax: +49 7182 93777 99
INDIA Tel: +91 80 4123 6041 Fax: +91 80 4123 6044 **IRISH REPUBLIC** Tel: 1800 409 011 **ITALY** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709
JAPAN Tel: +81 3 3661 7205 Fax: +81 3 3661 7206 **MIDDLE EAST** Tel: +9 71 4 360 0849 Fax: +9 71 4 368 8110 **SPAIN** Tel: +34 902 013 346 Fax: +34 902 013 379
SWITZERLAND Tel: 0800 83 47 63 Fax: 0800 83 47 54 **UK** (EU Authorised Representative) Tel: +44 1628 626 136 Fax: +44 1628 626 146 **Nthn Ireland** Tel: 0800 132 189
USA Tel: 1800 446 3908 or +1 949 470 3900 Fax: +1 949 470 3933